Kneehab<sup>XP</sup>

Distributed by:  
Neurotech Inc.  
A wholly-owned subsidiary of Theragen™ Inc.  
11220 Aztec Loop, Suite 101  
Montgomery, AL 36109  
USA  
Tel: 888-940-1997  
[www.neurotechinc.com](http://www.neurotechinc.com)

Designed by & Manufactured for:  
Theragen™ Inc.  
11220 Aztec Loop, Suite 101  
Montgomery, AL 36109  
USA

Copyright © 2012 Theragen™ Inc. All rights reserved.  
Caution: In the United States of America federal law restricts this device to sale or use by or on the  
order of a physician or other practitioner licensed by the laws of the state in which he/she practices.  
neurotech® and Kneehab™ XP are registered trademarks of Theragen™ Inc.

neurotech.  
embrace life.

Kneehab<sup>XP</sup>  
Quadriceps Therapy System

CONTROLLER

TYPE 412

## INSTRUCTIONS FOR USE

Part No. 3401-0305  
Rev. 5  
Issued: 5/19

Kneehab<sup>XP</sup>

## Kneehab+ Controller Instructional Photos


Figure A

Figure B

Figure C

Figure D

Figure E

Figure F

Figure G

## ENGLISH – INTRODUCTION & CONTACT DETAILS

### Introduction

Dear Customer

Thank you for choosing Kneehab® XP Neuromuscular Electrical Stimulator (NMES) Controller. Neurotech® has developed Kneehab XP to provide a uniquely convenient way to prevent or retard disuse atrophy in the quadriceps (thigh muscle weakness). Kneehab XP is designed to deliver strong, comfortable quadriceps contractions thereby ensuring a faster rehabilitation period. It utilizes Multipath™ stimulation technology, which is a highly advanced form of muscle stimulation for injury rehabilitation.

To see the benefits of Kneehab XP and Multipath technology in more detail, please see the Frequently Asked Questions on page 7 for a full explanation.

The Kneehab XP also includes Transcutaneous Electrical Nerve Stimulation (TENS). TENS is a clinically proven, safe and effective means of providing a drug-free method of pain relief. Chronic pain sufferers who experience sudden breakthrough levels of knee pain, can rely on Kneehab XP for immediate and long-term pain relief when used with drug and non-drug treatment regimens.

Your Kneehab XP pack contains your Kneehab XP Controller, battery charger and a user manual.

If you have further questions regarding Kneehab XP, please contact your prescribing physician, the distributor who provided you with Kneehab XP or Neurotech at the address/phone number below.

Yours sincerely,  
The Neurotech Team.

### If you have questions or require further information please contact:

**Neurotech NA Inc.**  
A wholly-owned subsidiary of Theragen Inc.  
11220 Assett Loop, Suite 101  
Manassas VA 20109  
USA  
Tel: 888-980-1197  
Email: patientcare@neurotech.us  
[www.kneehab.us](http://www.kneehab.us)

# TABLE OF CONTENTS

---

<b>1.0</b>	<b>Safety Information</b>	<b>4</b>
<b>2.0</b>	<b>TENS and NMES Information</b>	<b>4</b>
2.1	Contraindications -TENS	5
2.2	Warnings -TENS	5
2.3	Precautions: - TENS	6
2.4	Adverse Reactions: TENS	6
2.5	Contraindications - NMES	7
2.6	Warnings - NMES	7
2.7	Precautions: - NMES	8
<b>3.0</b>	<b>Frequently Asked Questions</b>	<b>9</b>
<b>4.0</b>	<b>Getting to Know Kneehab XP</b>	<b>11</b>
4.1	Kneehab XP Pack Contents	11
4.2	Description of Controls	12
4.3	Description of Unit Display	13
<b>5.0</b>	<b>Operating Kneehab XP</b>	<b>13</b>
5.1	Problem Solving Guide	15
<b>6.0</b>	<b>Program Information</b>	<b>16</b>
<b>7.0</b>	<b>System Maintenance</b>	<b>17</b>
7.1	Accessories	18
<b>8.0</b>	<b>Warranty</b>	<b>18</b>
<b>9.0</b>	<b>Technical Specifications</b>	<b>19</b>

## 1.0 Safety Information

### Validity

The information and technical data contained in this document relates to the Kneehab XP Controller provided with this User Manual. Each Kneehab XP Controller is attributed a serial number which is located on the back of the Controller. The information and technical data disclosed in this document are proprietary to Neurotech and may only be used and disseminated for the purposes and to the extent specifically authorized in writing by the company.

### Disclaimers

All items of equipment manufactured and sold by Neurotech, a Division of Theragen Inc., are rigorously checked and tested prior to shipment. However, the use of this equipment is outside of the control of the company. Neurotech only accepts responsibility for the safety, reliability and performance of the equipment when it is operated in accordance with the instructions herein and within the given specifications. Therefore, the user must bear full responsibility for any actions arising out of the use or misuse of this equipment. Any modifications, repairs or servicing must be undertaken by authorized Neurotech personnel.

Kneehab XP is designed by and manufactured for Neurotech. Copyright © 2012 by Theragen Inc., all rights reserved. neurotech® is a registered trademark of Theragen Inc.

The sale and/or operation of this equipment is subject to legislation in a number of localities. Compliance with this legislation rests with the dealer, or user of the equipment as appropriate.

Clinical trials of Kneehab XP in the post-surgical setting are currently limited to ACL surgery.

### Prescription Required:

Caution: In the United States of America federal law restricts the device to sale or use by, or on the order of a physician or other practitioner licensed by the laws of the state in which he/she practices.

## 2.0 NMES and TENS Information

### Kneehab XP - Transcutaneous Electrical Nerve Stimulation (TENS) Capabilities:

#### Intended Use:

Kneehab XP delivers stimulation based on the principles of transcutaneous electrical nerve stimulation (TENS) as described in the following. Short electrical pulses are sent via self-adhesive conductive gel pads to the surface of the skin.

Powered muscle stimulators should only be used under medical supervision for adjunctive therapy for the indications listed below.

**TENS Indications:**

- Symptomatic relief and management of chronic, intractable pain;
- Adjunctive treatment in the management of acute, post-surgical or post-traumatic pain;
- Provide symptomatic relief and management of intractable pain and relief of pain associated with arthritis; and
- Provide an adjunctive therapy in reducing the level of pain and symptoms associated with osteoarthritis of the knee.

**Transcutaneous Electrical Nerve Stimulation (TENS)**

TENS is a pain therapy based on the application of electrical stimuli to the skin via stimulation of the nerve fibers. This comfortable, tingling or pulsing sensation is designed to suppress pain signals. It can also provide a carry over effect by encouraging the release of endorphins, the body's natural pain-relieving chemicals. TENS can be used to complement existing drug therapies for pain relief.

**2.1 Contraindications: TENS**

- Do not use this device on patients whose pain syndromes are undiagnosed; and
- Patients with electronic implants (e.g. cardiac pacemaker or defibrillator - as your product may interfere with the proper functioning of the implanted stimulator) or if you suffer from any other heart problem.

**2.2 Warnings - TENS**

- Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use;
- Do not apply stimulation when the patient is in the bath or shower;
- Do not apply stimulation while the patient is sleeping, unless the Controller keylock feature has been engaged;
- Do not apply stimulation while the patient is driving, operating machinery, or during any activity in which electrical stimulation can put the patient at risk of injury;
- Consult with the patient's physician before using this device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals;
- Apply stimulation only to normal, intact, clean, healthy skin;
- The long-term effects of chronic electrical stimulation are unknown;
- Stimulation should not be applied over the carotid sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex;
- Stimulation should not be applied over the neck or mouth. Severe spasm of the laryngeal and pharyngeal muscles may occur and the contractions may be strong enough to close the airway or cause difficulty in breathing;

- Stimulation should not be applied transthoracically in that the introduction of electrical current into the heart may cause cardiac arrhythmias;
- Stimulation should not be applied transcerebrally;
- Stimulation should not be applied over swollen, infected, or inflamed areas or skin eruptions, e.g., phlebitis, thrombophlebitis, varicose veins, etc; and
- Stimulation should not be applied over, or in proximity to, cancerous lesions.

## 2.3 Precautions: TENS

- TENS is not effective for pain of central origin, including headache;
- TENS is not a substitute for pain medications and other pain management therapies;
- TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism;
- Effectiveness is highly dependent upon patient selection by a practitioner qualified in the management of pain patients;
- Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head, and electrodes should not be placed on opposite sides of the head;
- The safety of electrical stimulation during pregnancy has not been established;
- Patients with suspected or diagnosed heart disease should follow precautions recommended by their physicians; and
- Patients with suspected or diagnosed epilepsy should follow precautions recommended by their physicians;
- Use caution when the patient has a tendency to bleed internally, such as following an injury or fracture;
- Use caution following recent surgical procedures when stimulation may disrupt the patient's healing process;
- Use caution if stimulation is applied over the menstruating or pregnant uterus;
- Use caution if stimulation is applied over areas of skin that lack normal sensation;
- Keep this device out of the reach of children;
- Use this device only with the leads, electrodes, and accessories recommended by the manufacturer; and
- Use this device only under the continued supervision of a licensed practitioner.

## 2.4 Adverse Reactions: TENS

- Patients may experience skin irritation and burns beneath the stimulation electrodes applied to the skin;
- Patients may experience headache and other painful sensations during or following the application of electrical stimulation near the eyes and to the head and face; and
- Patients should stop using the device and should consult with their physicians if they experience adverse reactions from the device.

## Neuromuscular Electrical Stimulation (NMES) Capabilities

### Intended Use:

Kneehab XP applies muscle and nerve stimulation using the principles of Neuromuscular Electrical Nerve Stimulation (NMES), which is defined below. It applies brief electrical pulses through skin surface adhesive conductive gel pads (available separately). Powered muscle stimulators should only be used under medical supervision for adjunctive therapy for the indications listed below.

### NMES Indications

- Maintain or increase the range of motion
- Prevention or retardation of disuse atrophy
- Re-educate muscles
- Early post-surgical quadriceps strengthening and improved post surgical knee stability secondary to quadriceps strengthening
- Relax muscle spasms
- Increase local blood circulation

### Neuromuscular Electrical Stimulation (NMES)

NMES may be defined as the application of electrical stimulation of the peripheral nervous system to contract a muscle either through the direct activation of the motor neurons in the mixed peripheral nerve, or indirectly through reflex recruitment.

## 2.5 Contraindications - NMES

- Powered muscle stimulators should not be used on patients with cardiac demand pacemakers/defibrillators.

## 2.6 Warnings - NMES

- The long-term effects of chronic electrical stimulation are unknown;
- Stimulation should not be applied over the carotid sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex;
- Stimulation should not be applied over the neck or mouth. Severe spasm of the laryngeal and pharyngeal muscles may occur and the contractions may be strong enough to close the airway or cause difficulty in breathing;
- Stimulation should not be applied transthoracically, in that the introduction of electrical current into the heart may cause cardiac arrhythmias;
- Stimulation should not be applied transcerebrally;
- Stimulation should not be applied over swollen, infected, or inflamed areas or skin eruptions, e.g. phlebitis, thrombophlebitis, varicose veins, etc; and
- Stimulation should not be applied over, or in proximity to, cancerous lesions.

## 2.7 Precautions - NMES

- If in doubt, always seek medical advice;
- Medical advice must be obtained before use on persons who are insulin-dependent diabetics or for persons who are under medical supervision for any cognitive dysfunction;
- Medical opinion must be obtained before persons with any serious illness or injury apply muscle stimulation;
- Powered muscle stimulators should be kept out of the reach of children;
- The safety of powered muscle stimulators for use during pregnancy has not been established;
- Caution should be used for patients with suspected or diagnosed heart problems;
- Caution should be used for patients with suspected or diagnosed epilepsy;
- Caution should be used in the presence of the following:
  - When there is a tendency to hemorrhage following acute trauma or fracture;
  - Following recent surgical procedures when muscle contraction may disrupt the healing process;
  - Over the menstruating or pregnant uterus;
  - Over areas of the skin which lack normal sensation.
- Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium gel;
- Stimulation settings should be based on the guidance of the prescribing practitioner;
- Powered muscle stimulators should be used only with the leads and conductive gel pads recommended for use by the manufacturer;
- Portable powered muscle stimulators should not be used while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk or injury;
- Operation in close proximity to shortwave or microwave therapy equipment may produce instability in the stimulator output;
- Electronic monitoring equipment (such as ECG monitors and ECG alarms) may not operate properly when the stimulator is in use.
- Although compliant with applicable EMC requirements, this device may still interfere with more sensitive equipment, please move away or switch off;
- Discontinue using the device following knee surgery if you experience any of the following adverse reactions: rotational trauma or hamstring strain during leisure time activities/sports, spasms of the quadriceps muscle, excessive post-operative swelling, redness and an evidence for haemorrhage; and
- In clinical tests of Kneehab XP post ACL surgery, patients with the following conditions were excluded: multi-ligament construction, additional meniscectomy, quadriceps tendon graft, cartilage damage greater than grade II, micro-fracture of the knee. Accordingly, safety and efficacy of the device with these patients has not been demonstrated.

**Do not use Kneehab XP if:**

- Positioning the conductive gel pads necessitates placement over areas at which drugs/ medicines are administered by injection (short term or long term) e.g. Hormone treatment.
- You are connected to high-frequency surgical equipment.

**Adverse Reactions**

- Skin irritation and burns beneath the conductive gel pads have been reported with the use of powered muscle stimulators. Contact your prescribing physician, or Neurotech if any irritation, skin reaction, hypersensitivity or other adverse reaction is experienced and discontinue using in the meantime. Note, however, that a slight reddening of the skin is quite normal under the conductive gel pads during and for a short time after treatment.
- Patients should stop using the device and should consult with their physicians if they experience adverse reactions from the device.

**Note:**

- Kneehab XP is designed for single patient use only.

**Important**

- The safety information provided in this User Manual must be followed.

## 3.0 Frequently Asked Questions

**Q: What is the Kneehab XP Controller?**

A: Kneehab XP Controller is a unique, innovative Neuromuscular Electrical Stimulator (NMES) product designed specifically for the prevention or retardation of disuse atrophy in the quadriceps (thigh muscle weakness). Kneehab XP can be used both pre- and post-operatively for the prevention or retardation of disuse atrophy in the quadriceps. Kneehab XP is clinically proven to provide early post-surgical quadriceps strengthening and improved post-surgical knee stability, secondary to quadriceps strengthening.

**Q: What is Multipath?**

A: Kneehab XP utilizes Multipath stimulation technology, which is a highly advanced form of NMES for injury rehabilitation. Multipath optimizes rehabilitation through strong, comfortable quadriceps contractions. EMS has been used for over 50 years to re-activate and strengthen weak muscles. Multipath is a major breakthrough in this area, utilizing sophisticated modern engineering and software technology to focus stimulation precisely in those parts of the muscle that optimize muscle contraction. Multipath stimulation enables Kneehab XP to deliver highly-focused, accurate and comfortable quadriceps contractions.

**Q: What should Kneehab XP feel like when operating?**

A: When you switch Kneehab XP on and gradually increase the intensity, you will feel a mild, comfortable “tingling” or “buzzing” sensation in your quadriceps muscles. As you increase the intensity, the muscles will begin to contract.

The channel on the inside of your thigh will increase the intensity in the inner part of the quadriceps (called the Vastus Medialis). For some conditions, your clinician may recommend you only use this channel. The channel on the outside of your thigh will increase the intensity in the rest of the quadriceps. Both channels are used together in general use.

Kneehab XP assists in the stabilization of the knee by contracting the inner quadriceps before the outer quadriceps.

**Q: How can Kneehab XP help rehabilitate the knee?**

A: The quadriceps muscle is the most important stabilizer of the knee. After knee injuries or knee surgery, pain and swelling often cause disuse, which leads to muscle atrophy. An essential component of rehabilitation after knee injury is strengthening of the quadriceps muscle. Kneehab XP works by re-educating and strengthening the quadriceps muscle through cycles of contraction and relaxation. Kneehab XP is a highly effective, yet simple and convenient way to rehabilitate quadriceps atrophy after injury or surgery.

**Q: How often should I use Kneehab XP?**

A: Your prescribing physician will advise you on the duration of your treatment program. The only treatment regimen that has been studied involves 12 weeks of therapy: 5 days per week, 3 times per day, 20 minutes each session.

**Q: What should I do if treatment becomes uncomfortable?**

A: Effective treatment should not cause undue discomfort. If treatment becomes painful, try lowering the intensity level slightly, then gradually build back up to previous level. If discomfort persists, Neurotech recommends that you contact your prescribing physician who will advise you on usage of Kneehab XP to minimize discomfort.

**Q: Can Kneehab XP be used in conjunction with other forms of therapy?**

A: Kneehab XP can be used as an adjunct to most forms of therapy for knee injuries. See page 5 of the User Manual for Contraindications (situations when Kneehab XP should not be used).

**Q: Can I move around while using Kneehab XP?**

A: If you need to move around during a treatment session, press the On/Off (Pause) button briefly, to pause the treatment. To complete the treatment session, just press the On/Off (Pause) button a second time.

**Q: Who should I contact if I have additional questions?**

A: Contact your prescribing physician or, the distributor that provided you with Kneehab XP or Neurotech if you have any questions regarding usage of Kneehab XP.

**Q: What should I do if I experience an adverse reaction?**

A: If any irritation, skin reaction, hypersensitivity or other adverse reaction is experienced, discontinue use and contact your prescribing physician. A slight reddening of the skin is normal under the conductive gel pads during and for a short time after treatment. This should disappear within approximately 30 minutes.

## 4.0 Getting to know Kneehab XP

Kneehab XP Controller is a portable, two-channel neuromuscular electrical stimulator (NMES), which operates using constant current pulses to stimulate the nerves in the quadriceps area of the body using Neurotech's Multipath technology. These pulses are designed to cause muscular contractions through the application of electrical stimulation to the peripheral nervous system. Kneehab XP is an easy to use solution, specifically designed for home use.

**Battery Information**

The unit is powered by a 3.6V NiMH rechargeable battery. A battery charger is included with the product.

### 4.1 Kneehab XP Pack Contents

**1. Kneehab XP Controller**

The unit generates signals which are sent, via the conductive gel pads, to your muscles, to make them contract. Your unit is rechargeable and takes approximately 3 hours to charge fully.

**2. Battery Charger:**

Kneehab XP is supplied with a charger. To charge the battery pack, connect the Controller to the charger and plug the charger into a main's supply socket.

**3. Kneehab XP User Manual:**

The Kneehab XP User Manual provides a guide to product use and safety information.

## 4.2 Description of Controller Controls

### 1. On/Off Button (○)

Press and hold this button for 3 seconds to turn your unit on or off. During a treatment session, you can pause the stimulation by briefly pressing this button. To un-pause, briefly press this button again.

### 2. Intensity Buttons (▼ / ▲)

Press and hold these buttons to increase and decrease the intensity of the stimulation. The right-hand buttons control the right-hand side of the garment. The left-hand buttons control the left-hand side of the garment.

### 3. Program Button (P)

Press and hold this button for 3 seconds before starting a treatment session to change the selected stimulation program. **Note:** See section 5.1 for program details.

Press and hold this button for 3 seconds during a treatment session to activate/ deactivate the intensity button lock function. You can use this function to prevent accidental changes of intensity during a treatment session. Use the lock function if you plan to sleep during treatment.


### 4. Information Button (i)

After treatment has started:

- Once treatment has started and you are using an NMES program, this button will no longer function as an information button but instead will operate as a trigger function.
- The trigger function is used for stimulating a muscle contraction on-demand during a therapy program. Activate the trigger function by pressing the information button (i) during an NMES program (programs #1-#6), it does not work during TENS (programs #7-#9). The muscle contraction is sustained for as long as the (i) button is held down. The muscle contraction stops when the button is released. The muscle will contract again if the (i) button is pressed. On the screen, the trigger function is signaled by 2 flashing side bars. The trigger function can be deactivated by pressing any intensity button, up or down. The KneeHab XP Controller returns to the preset program once the trigger function is deactivated.

Before starting the programs, press this button briefly at any time to toggle through the following displays:

- Current Intensity level for left and right sides
- Total Treatment Time in hours and minutes

- Average Intensity level over the last 4 sessions for left and right sides
- Return to original display

Before treatment has started, you may press and hold this button for 3 seconds to activate/ deactivate the mute function. The ( i ) button can also be used to reset the Controller's timer. To reset timer, press the information key ( i ) and the right down amplitude key together for 2 seconds.

### **4.3 Description of Controller Display**

The following icons appear on the Controller display at various times throughout the treatment.

Toning intensity from the left-hand side of the garment

Toning intensity from the right-hand side of the garment

Counts down the time left in the current session in minutes and seconds.

Indicates that the session has been paused.

Indicates that the mute function is activated

Intensity button lock is active

This warning symbol indicates a poor contact between the unit and garment or between the conductive gel pads and the skin (see page 13 for more details).

Shows the battery power remaining.

Shows the program currently running.

This symbol appears around the edge of the display during each contraction phase.

## **5.0 Operating Kneehab XP Controller**

### **Step 1 - Charge the battery**

Connect the Controller to the charger (Figure A) and plug the charger into an electrical outlet. To remove the Controller from the charger press the buttons on each side of the unit connector and pull apart. You will know the battery pack is on charge when the unit display shows a battery symbol with 3 bars that turn on and off. It will take approximately 2-3 hours to fully charge. The battery pack is fully charged when the battery symbol and its 3 bars remain on.

**Notes:**

- Do not leave your Controller connected to the charger when the battery is fully charged.
- Do not leave your charger plugged into electrical outlet when not in use.

**Step 2 – Connect the Controller and switch it on**

Connect the Controller to an appropriate conductive garment by plugging it into the garment socket (Figure B). Switch on the Controller by pressing the On/ Off button for 3 seconds (Figure C). The Controller display will show the program number selected, the duration of the program and the battery status.

The Kneehab XP is designed for use on the Right leg (R) or Left Leg (L). The garment is marked on the outside, top with the letter R or L. Before starting treatment, confirm the letter printed on your garment is correct per the prescription from your health care professional. You must use correct garment on corresponding leg. When connected to a Left garment the Controller will display 'L' for 3 seconds at start up. When connected to a Right garment the Controller will display 'r'. Please check that your Controller displays the correct code for the garment to which it's connected. If the Controller does not display the expected code, please contact your distributor.

**Step 3 – Begin treatment**

Begin increasing the intensity (Figure D). This starts the treatment session and the timer will begin to count down. As you increase the left or right intensity button you will feel the stimulation increasing on the corresponding side of your leg. You will also see the intensity level for that side increasing on the Controller's display. Increase the intensity on both sides of your leg until you feel a firm, comfortable contraction of the quadriceps. The intensity level ranges from 0 - 99.

**Important Notes:**

- Typically, you should increase the intensity to the inside of your leg first and to optimize contractions, the intensity level for the inner thigh should be set at a higher level than the outer thigh. As you become more experienced using Kneehab XP, you will quickly find the levels which suit you best.
- The treatment is delivered in contraction and relaxation cycles. See your Program Information Card for more details.
- Changing the intensity at any time will trigger the start of a contraction cycle.
- Caution is required when using Kneehab XP because the strong quadriceps contraction it elicits may cause involuntary lower leg movement. If not otherwise directed by your prescribing physician, we recommend using Kneehab XP while sitting with your leg bent at an angle between 0-60 degrees (fig d) and your foot secured firmly to prevent lower leg movement.

## Step 4 - Finishing Treatment

The treatment session is completed when the counter reaches zero. The Controller will beep to indicate this. Switch the Controller off by pressing the On/ Off button for 3 seconds. Should you forget to do this, the Controller will switch off automatically after 20 seconds anyway.

Disconnect the Controller from the garment by pressing the buttons on either side of the connector and pulling the two sections apart (Figure E).

## 5.1 Problem Solving Guide

Problem	Possible cause	Solution
The display doesn't come on and there's no signal from the Controller	<ul style="list-style-type: none"> <li>• Controller not switched on</li> <li>• Battery pack is fully discharged</li> </ul>	<ul style="list-style-type: none"> <li>• Press and hold the On/Off button for 3 seconds to switch the Controller on</li> <li>• Recharge the battery pack</li> </ul>
The battery symbol is flashing/ there's ineffective stimulation	<ul style="list-style-type: none"> <li>• The battery charge is low</li> </ul>	<ul style="list-style-type: none"> <li>• Recharge the battery pack</li> </ul>
The contractions are very weak, even when the intensity is high	<ul style="list-style-type: none"> <li>• The battery charge is low</li> </ul>	<ul style="list-style-type: none"> <li>• Recharge the battery pack</li> </ul>
The warning symbol appears on the display 	<ul style="list-style-type: none"> <li>• Controller is not properly attached to the garment</li> </ul>	<ul style="list-style-type: none"> <li>• Ensure the Controller and garment are properly connected</li> </ul>
ERR has appeared on the display	<ul style="list-style-type: none"> <li>• Error message, indicating a problem with the Controller</li> </ul>	<ul style="list-style-type: none"> <li>• Switch the Controller off and on again. It should now work without problems. If the error message persists, call your distributor</li> </ul>

## 6.0 Program Information

### Neuromuscular Electrical Stimulation (NMES) Programs on Kneehab XP:

Program Number (Note 2) (on Screen)	Duration (Minutes)	Frequency/ Rate (Hz or pps)	Pulse Width (usec)	Ramp-Up time (seconds)	Contraction Time (seconds)	Ramp Down Time (seconds)	Relaxation Time (seconds)	Additional Function	Indication(s) (See Note 1)
P1	20	50Hz	300 $\mu$ s– 400 $\mu$ s	1	5	0.5	10	Trigger	4
P2	20	50Hz	300 $\mu$ s– 400 $\mu$ s	1	10	0.5	10	Trigger	1,2,3,5,6
P3	20	50Hz	300 $\mu$ s– 400 $\mu$ s	1	10	0.5	20	Trigger	1,2,3,5,6
P4	20	50Hz	300 $\mu$ s– 400 $\mu$ s	1	10	0.5	30	Trigger	1,2,3,5,6
P5	20	35Hz	300 $\mu$ s– 400 $\mu$ s	1	5	0.5	5	Trigger	1,2,3,5,6
P6	20	70Hz	300 $\mu$ s– 400 $\mu$ s	1	10	0.5	50	Trigger	1,2,3,5,6

**Note 1:**

**Indications are numbered as follows:**

1. Maintain or increase the range of motion.
2. Prevention or retardation of disuse atrophy.
3. Re-educate muscles.
4. Early post-surgical quadriceps strengthening and improved post surgical knee stability secondary to quadriceps strengthening.
5. Relax muscle spasms.
6. Increase local blood circulation.

**Note 2:**

These programs use Multipath technology.

### Transcutaneous Electrical Nerve Stimulation (TENS) Programs on Kneehab XP:

Program Number (on Screen)	Duration (Minutes)	Frequency/ Rate (Hz or pps)	Pulse Width (usec)	Ramp-Up Time (seconds)	Contraction Time (seconds)	Ramp Down Time (seconds)	Relaxation Time (seconds)	Additional Function	Indication(s) (See Note 3 on page 17)
P7	30	99Hz	200 $\mu$ s	N/A	Continuous	N/A	N/A	No-Trigger	7 to 10
P8	30	4Hz	300 $\mu$ s	N/A	Continuous	N/A	N/A	No-Trigger	7 to 9
P9	30	125Hz	175 $\mu$ s	N/A	Continuous	N/A	N/A	No-Trigger	7 to 9

**Note 3:****Indications are numbered as follows:**

7. Symptomatic relief and management of chronic, intractable pain.
8. Adjunctive treatment in the management of acute, post-surgical or post-traumatic pain.
9. Provide symptomatic relief and management of intractable pain and relief of pain associated with arthritis.
10. Provide an adjunctive therapy in reducing the level of pain and symptoms associated with osteoarthritis of the knee.

## 7.0 System Maintenance

### Battery Power / Replacing the Battery

The battery icon on the display will flash when the battery is low and needs to be recharged. After a period of time, you may find that the battery will only power your Controller for a limited number of sessions, requiring you to charge it more frequently. This is because rechargeable batteries have a limited number of charge cycles and may eventually need to be replaced.

Should you need to replace the rechargeable battery pack, remove the small rubber plug from the back cover of your Controller (Figure F), unscrew the back cover and remove it. Replace the existing battery pack with a new battery pack (Figure G). Replace the battery cover and rubber plug. You can purchase a new battery pack by contacting the distributor who provided you with this product.

### Repair, Service & Modification

Access to the interior is not required for maintenance purposes. If your Controller is damaged, you should not use it, but should return it to Neurotech or your local distributor for replacement or repair. Repairs, service and modifications may not be carried out by anyone other than qualified service personnel authorized by Neurotech. Neurotech will not accept any responsibility if these guidelines and instructions are not followed. Only parts specified by Theragen (available from Neurotech or Theragen) should be used with the Kneehab XP.

### Cleaning Instructions

You should never let your Kneehab XP Controller get wet, but it may be cleaned using a soft cloth, lightly dampened in mild soapy water.

## 7.1 Accessories

**Kneehab XP Conductive Garment:** Right Leg (1274-3311)  
Left Leg (1274-3321)

**Kneehab XP Battery Pack:** (3.6V NiMH)

**Important:** Under no circumstance should anything other than the correct type of battery pack (rechargeable, 3.6V NiMH) be used with your Controller. These can be purchased from the Neurotech office.

### Disposal of battery packs

Used battery packs must never be disposed of in a fire, but in accordance with your state's laws governing the disposal of such items. Return all other contents to Neurotech, Theragen or the distributor who provided you with Kneehab XP.

### Battery Charger (US/Japan - 2504-0322)

Charger complies to EN 60950, Input: nominal voltage 100-240V, frequency 50-60Hz, current 125mA. Do not use any other chargers or any other power supplies.

### Kneehab XP Electrodes or Conductive Gel Pads:

Type: 613/ 614/ 615/ 616  
Kneehab XP uses electrodes or conductive gel pads which have been specifically designed for the Kneehab XP. Replacement adhesive electrodes or conductive gel pads are available. To order replacement electrodes or conductive gel pads, please contact the distributor that provided you with Kneehab XP. Only conductive gel pads specified by Neurotech. for use with Kneehab XP may be used. Using other electrodes or conductive gel pads may degrade performance levels. All electrodes or conductive gel pads are provided with their own instructions where necessary, including cleaning, storage and care instructions. For hygiene reasons, the garment and conductive gel pads are for single person use only. With proper care, as per this manual and typical use, the gel pads have an estimated usage of approximately one month.

Note: Conductive gel pads wear out over time: If they are dirty or no longer adhere properly, they need to be replaced. Replace the leads or garment if the sheathing is damaged and exposes the copper wire.

## 8.0 Warranty

**neurotech®** offers a warranty of 3 years for the Kneehab XP Controller\*. Should your Kneehab XP develop a fault within the specified time, Neurotech will replace or repair your Kneehab XP free of charge, provided the unit:

- has been used for its intended purpose and in the manner described in this instruction manual.
- has not been connected to an unsuitable power source.

- has not been subjected to misuse or neglect.
- has not been modified or repaired by anyone other than an approved **neurotech®** distributor or agent.

This warranty complements existing national guarantee obligations and does not affect your statutory rights as a consumer.

\* excludes battery pack.

## 9.0 Technical Information

### General Specifications:

**Product Type:** 421

**Classification:** Internally powered equipment, Class II charger, Type BF applied parts.

**Intended use:** Electrical Muscle Stimulator

**Waveform:** Symmetrical bi-phasic square waveform when measured into a standard bodyload.

### Environmental Specifications:

**Operating Range:** Temperature: 10 to 35°C / 55 to 95°F

Humidity: 20 to 65% RH

**Transport & Storage Range:** Temperature: 0 to 55°C / 32 to 131°F

Humidity: 10 to 90% RH

### Electrical Specifications:

#### Rated Outputs - Voltage / Currents

Parameter	500 Ω	1K Ω	1.5K Ω
Output RMSV	11V	22V	26.5V
Output RMSA	22mA	22mA	17.5mA
Output Frequency	4 - 125Hz	4 - 125Hz	4 - 125Hz
DC Component: Approx	0 C	0 C	0 C
Positive Pulse Width	100 - 400 µS	100 - 400 µS	100 - 400 µS
Negative Pulse Width	100 - 400 µS	100 - 400 µS	100 - 400 µS
Interphase Interval	40 - 100 µS	40 - 100 µS	40 - 100 µS

## Description of your Controller's symbols:

There are a number of technical markings on your Controller. These can be explained as follows:

The unit requires 1 x 3.6 volt (NiMH) DC battery pack DC is indicated by the symbol: ===

Output (RMSA) Stands for the effective current output, which is the root mean square current measured at a specified resistance.

Output (RMSV) Stands for the effective voltage output, which is the root mean square voltage measured at a specified resistance.

Power (P): Maximum power output measured in Watts (W) into a  $500\ \Omega$  load.

Frequency (F): Number of pulses output by the unit per second, measured in Hertz (Hz).

 This symbol means "Attention, consult the accompanying documents".

 This symbol means type BF equipment.

 This symbol means "Keep Dry"

 This symbol means "Read Instructions"

 This symbol means "Recyclable Materials"

SN stands for "serial number".

On the Kneehab XP Controller is the serial number specific to this unit. The letter preceding the number indicates the year of manufacture, where "L" denotes 2006, "M" denotes 2007 etc.

Conforms to UL STD 60601-1.  
Certified to CAN/CSA.  
STD C22.2 NO 601.1

### Disposal of device

At the end of the product lifecycle, do not throw out this product with the normal household garbage, but bring it to a collection point for the recycling of electronic equipment.

Some product materials can be re-used if you bring them to a recycling point. By re-using some parts or raw materials from used products you can make an important contribution to the protection of the environment. Please contact your local authorities if you need more information about collection points in your area.

Waste Electrical and Electronic Equipment can have potentially harmful effects on the environment. Incorrect disposal can cause harmful toxins to build up in the air, water and soil and can be harmful to human health.

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Kneehab XP est conçu et fabriqué par Neurotech. Copyright ©2012 par Theragen Inc., tous droits réservés neurotech® est une marque de commerce déposée de Theragen Inc.
Kneehab XP ha sido diseñado y manufacturado por Neurotech. Copyright© 2012 de Theragen Inc., derechos reservados. neurotech® es una marca registrada de Theragen Inc.